


**INFORMATIVE NOTE**
**Name: REFLEXX 40**
**Date: 04/01/2024 Rev: 018**
**MANUFACTURER**

 REFLEXX SpA  
 Via Passeri, 2 46019 VIADANA (MN)  
 Tel.: + 39 0375 758891 Fax + 39 0375 464504 [info@reflexx.com](mailto:info@reflexx.com)  
[www.reflexx.com](http://www.reflexx.com)

**Imballo - Pack**

100		10	
Taglia Size	Articolo Article	Taglia Size	Cartone Carton
R40 XS	8032891631324	XS	8032891636329
R40 S	8032891631331	S	8032891636336
R40 M	8032891631348	M	8032891636343
R40 L	8032891631355	L	8032891636350
R40 XL	8032891631362	XL	8032891636367

**DESCRIPTION:** ■ disposable examination glove powdered, not sterile ■ **Material:** natural latex ■ **Shape:** ambidextrous; five fingers, beaded cuff ■ **Colour:** white ■ **Feature:** smooth surface ■ **Packing:** 100 pieces/box; 10 boxes/carton ■ **ISO 2859:** AQL 1.5 (Level G1) freedom from holes; AQL 2.5 (level G1) major defects; AQL 4 (level G1) minor

**CONFORMITY AND PERFORMANCES**
**Personal Protective Equipment cat III**

0465 ANCI\_CIMAC, Via Aguzzafame, 60/b 27029 Vigevano (PV)

Reg Eu 2016/425 - EN ISO 374-1:2016 TIPO C - EN ISO 374-5:2016 - EN 374-2:2014 - EN 420:2003 + A1 2009

**Medical Device**

Reg EU 2017/745 - EN 455 1,2,3 &amp; 4

(CND T010201 NON SURGICAL LATEX GLOVES)



Reg CE 1935/2004 – Reg EU 10/2011 - DM 21/03/73. To such use carefully read the declaration of suitability for food contact (DICA).

DIMENSIONS		
Cod./Size	Length min* (mm)	Width. (±5 mm)
R40/ XS (6)		
R40/ S (7)	240	85
R40/ M (8)	240	95
R40/ L (9)	240	106
R40/ XL (10)	240	116

\* Gloves for special applications 5.1.3

<b>Thickness mm</b>	Palm (+/- 0.02)	0.07
---------------------	-----------------	------

<b>EN 455-2:2009/A1</b>	Tensile strength	≥6 N
-------------------------	------------------	------

<b>EN ISO 374-5:2016</b>	Protection from microorganisms (bacteria and fungi)	Virus protection ISO 16604:2004
	pass	pass

**REG EU 2016/425**

EN ISO 374-1:2016 + A1:2018/Type C


**EN ISO 374-5: 2016**


EN ISO 374-1:2016				
Chemical compound:	COD	EN 16523-1:2015 TYPE C		EN 374-4:2013
		Performance level	Breakthrough time	% degradation
Sodium hydroxide 40%	K	3	>60'	5%

RISKS FROM WHICH THE DPI IS DESIGNED TO PROTECT:

COD	Protection
K Sodium hydroxide 40%	CAS:1310-73-2 Protection from use of inorganic bases

**RECOMENDATIONS OF USE** Multi-purpose glove to protect the hand and/or the handled product, can also be used for general uses ■ This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals ■ The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm – where the cuff is tested also) and

REFLEXX S.P.A. UNIPERSONALE VIA PASSERI, 2 - 46019 VIADANA (MN) | TEL. + 39 0375 758891

 Cap.Soc. 1.200.000 € (i.v.) | PIVA 02085450209 | M5UXCR1 | R.E.A. 223166 | [info@reflexx.com](mailto:info@reflexx.com) [www.reflexx.com](http://www.reflexx.com)


5010014617

relates only to the chemical tested. It can be different if the chemical is used in a mixture ■ It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation ■ When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves ■ Before usage, inspect the gloves for any defect or imperfections ■ The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen ■

**WARNINGS** Before using the product, try it in the real condition of use ■ The glove is a device subject to wear, so frequent replacement is recommended and in all cases where imperfections are observed ■ Wear gloves with dry, clean hands ■ The product contains natural latex which can cause allergic reactions of type I and IV hypersensitivity. In case of reaction, discontinue use and consult a doctor ■ The user must inform the manufacturer and the relative authority of any incidents occurring during the use of the protective equipment

**STORAGE AND DISPOSAL** Store the gloves in their original packaging in a dry place. ■ Do not expose to direct sunlight and heat sources. ■ The expiry date is 3 years from the date of manufacture indicated on the sales package, and is valid for the product properly stored ■ Gloves must be disposed of according to local regulations ■ Declaration of Conformity on [www.reflexx.com](http://www.reflexx.com)