


INFORMATIVE NOTE
Name: REFLEXX 76 soft
Date: 04/09/2024 Rev: 013
MANUFACTURER

 REFLEXX SpA
 Via Passeri, 2 46019 VIADANA (MN)
 Tel.: + 39 0375 758891 Fax + 39 0375 464504 info@reflexx.com
www.reflexx.com

Imballo - Pack

100		10	
Taglia Size	Articolo Article	Taglia Size	Cartone Carton
R76 XS	8032891631874	XS	8032891636879
R76 S	8032891631881	S	8032891636886
R76 M	8032891631898	M	8032891636893
R76 L	8032891631904	L	8032891636909
R76 XL	8032891631911	XL	8032891636916

DESCRIPTION: ■ disposable examination glove powder free, not sterile ■ **Material:** nitrile ■ **Shape:** ambidextrous; five fingers, beaded cuff ■ **Colour:** blue ■ **Feature:** internally chlorinated; textured surface and textured fingertips ■ composition soft and sensitive ■ **Packing:** 100 pieces/box; 10 boxes/carton ■ **ISO 2859:** AQL 1.5 (level G1) freedom from holes; AQL 2.5 (level G1) major defects; AQL 4 (level G1) minor defects

CONFORMITY AND PERFORMANCES
Personal Protective Equipment cat III

0465 ANCI_CIMAC, Via Aguzzafame, 60/b 27029 Vigevano (PV)

Reg EU 2016/425 - EN ISO 374-1:2016 TYPE C - EN ISO 374-5:2016 - EN 374-2:2014 - EN 420:2003 + A1 2009

Medical Device

Reg EU 2017/745 – EN 455 1,2,3 & 4

(CND T01020204 NON SURGICAL NITRILE GLOVES)



Reg CE 1935/2004 – Reg EU 10/2011 - DM 21/03/73. To such use carefully read the declaration of suitability for food contact (DICA).


DIMENSIONS		
Cod./Size	Length min(mm)*	Width mm (±5)
R76/ XS (6)	240	75
R76/ S (7)	240	85
R76/ M (8)	240	95
R76/ L (9)	240	105
R76/ XL (10)	240	115

* Gloves for special applications 5.1.3

Thickness mm	Palm (+/- 0.02)	0.07
---------------------	-----------------	------

EN 455-2:2009/A1	Tensile strength	≥6 N
-------------------------	------------------	------

EN ISO 374-5:2016	Protection from microorganisms (bacteria and fungi)	Virus protection ISO 16604:2004
	pass	pass

REG EU 2016/425


EN ISO 374-1:2016				
Chemical compound:	COD	EN 16523-1:2015 TYPE C		EN 374-4:2013
		Performance level	Breakthrough time	% degradation
Sodium hydroxide 40%	K	3	>60'	7,1%

RISKS FROM WHICH THE DPI IS DESIGNED TO PROTECT:

COD	Protection
K Sodium hydroxide 40%	CAS:1310-73-2 Protection from use of inorganic bases

RESULTS OF TESTS FOR ADDITIONAL CHEMICAL COMPOUNDS ACCORDING TO EN ISO 374-1: 2016 THE FOLLOWING RESULTS DO NOT GIVE PRESUMPTION OF CONFORMITY:

EN ISO 374-1:2016				
Chemical compound :	COD	EN 16523-1:2015		EN 374-4:2013
		Performance level	Breakthrough time	% degradation
Formaldehyde 37%	T	3	>60'	16,02%
Hydrogen peroxide 30%	P	1	>10'	-3,89%

RECOMENDATIONS OF USE Multi-purpose glove to protect the hand and/or the handled product, can also be used for general uses ■ The glove is suitable for contact with pure chemical products for a period of less than the permeation time indicated in the chemical resistance table ■ Results obtained under laboratory conditions may not reflect actual working time

WARNINGS Before using the product, try it in the real condition of use ■ The glove is a device subject to wear, so frequent replacement is recommended and in all cases where imperfections are observed ■ Wear gloves with dry, clean hands ■ In case of transient skin reactions and hypersensitivity reactions, discontinue use and consult a doctor ■ The user must inform the manufacturer and the relative authority of any incidents occurring during the use of the protective equipment

STORAGE AND DISPOSAL Store at room temperature ■ Store the gloves in their original packaging in a dry place. ■ Do not expose to direct sunlight and heat sources. ■ The expiry date is 3 years from the date of manufacture indicated on the sales package, and is valid for the product properly stored ■ Gloves must be disposed of according to local regulations ■ Declaration of Conformity on www.reflexx.com