


**INFORMATIVE NOTE**
**Name: REFLEXX 78**
**Date: 09/06/2023 Rev: 012**
**MANUFACTURER**

 REFLEXX SpA  
 Via Passeri, 2 46019 VIADANA (MN)  
 Tel.: + 39 0375 758891 Fax + 39 0375 464504  
[info@reflexx.com](mailto:info@reflexx.com)  
[www.reflexx.com](http://www.reflexx.com)


Cod./Size	Pack	
	Box	Carton
R78/XS	8032891633595	8032891638590
R78/S	8032891631928	8032891636923
R78/M	8032891631935	8032891636930
R78/L	8032891631942	8032891636947
R78/XL	8032891631959	8032891636954

**DESCRIPTION:** ■ disposable glove powder free, not sterile ■ **Material:** nitrile ■ **Shape:** ambidextrous; five fingers, beaded cuff ■ **Colour:** black ■ **Feature:** internally chlorinated; textured surface ■ **Packing:** 100 pieces/box; 10 boxes/carton ■ **ISO 2859:** AQL 1.5 (level G1) freedom from holes; AQL 2.5 (level G1) major defects; AQL 4 (level G1) minor defects

**CONFORMITY AND PERFORMANCES**
**Personal Protective Equipment cat III**

 0465 ANCI\_CIMAC, Via Aguzzafame, 60/b 27029 Vigevano (PV)  
 Reg EU 2016/425 - EN ISO 374-1:2016 TYPE C - EN ISO 374-5:2016 VIRUS (ISO 16604:2004). - EN 374-2:2014 - EN 420:2003 + A1 2009

**Medical Device**

 Reg EU 2017/745  
 (CND T01020204 NON SURGICAL NITRILE GLOVES)


Reg CE 1935/2004 – Reg EU 10/2011 - DM 21/03/73. To such use carefully read the declaration of suitability for food contact (DICA).

DIMENSIONS		
Cod./Size	Length min(mm)	Width (±5 mm)
R78/ XS (6)	240	75
R78/ S (7)	240	85
R78/ M (8)	240	95
R78/ L (9)	240	105
R78/ XL (10)	240	115

<b>Thickness mm</b>	Palm (+/- 0.02)	0.08
---------------------	-----------------	------

<b>EN 455-2:2009/A1</b>	Tensile strength	≥6 N
-------------------------	------------------	------

<b>EN ISO 374-5:2016</b>	Protection from microorganisms (bacteria and fungi)	Virus protection ISO 16604:2004
	pass	pass

**REG EU 2016/425**


EN ISO 374-1:2016				
Chemical compound :	COD	EN 16523-1:2015 TYPE C		EN 374-4:2013
		Performance level	Breakthrough time	% degradation
Sodium hydroxide 40%	K	3	>60'	7,1%

**RISKS FROM WHICH THE DPI IS DESIGNED TO PROTECT:**

Chemical compound :	Protection
K Sodium hydroxide 40 %	CAS:1310-73-2 Protection from use of inorganic bases

 RESULTS OF TESTS FOR ADDITIONAL CHEMICAL COMPOUNDS ACCORDING TO EN ISO 374-1: 2016  
 THE FOLLOWING RESULTS DO NOT GIVE PRESUMPTION OF CONFORMITY:

EN ISO 374-1:2016				
Chemical compound :	COD	EN 16523-1:2015		EN 374-4:2013
		Performance level	Breakthrough time	% degradation
Formaldehyde 37%	T	2	>30'	50,15%
Hydrogen peroxide 30%	P	1	>10'	15,19%

REFLEXX S.P.A. UNIPERSONALE VIA PASSERI, 2 - 46019 VIADANA (MN) | TEL. + 39 0375 758891

 Cap.Soc. 1.200.000 € (i.v.) | PIVA 02085450209 | M5UXCR1 | R.E.A. 223166 | [info@reflexx.com](mailto:info@reflexx.com) [www.reflexx.com](http://www.reflexx.com)


5010014617

**RECOMENDATIONS OF USE** Multi-purpose glove to protect the hand and/or the handled product, can also be used for general uses ■ The glove is suitable for contact with pure chemical products for a period of less than the permeation time indicated in the chemical resistance table ■ Results obtained under laboratory conditions may not reflect actual working time **WARNINGS** Before using the product, try it in the real condition of use ■ The glove is a device subject to wear, so frequent replacement is recommended and in all cases where imperfections are observed ■ Wear gloves with dry, clean hands ■ In rare cases where transient hypersensitive skin reactions may occur, discontinue use and consult a doctor ■ The user must inform the manufacturer and the relative authority of any incidents occurring during the use of the protective equipment **STORAGE AND DISPOSAL** Store the gloves in their original packaging in a dry place. ■ Do not expose to direct sunlight and heat sources. ■ The expiry date is 3 years from the date of manufacture indicated on the sales package, and is valid for the product properly stored ■ Gloves must be disposed of according to local regulations ■ Product Declaration of Conformity on [www.reflexx.com](http://www.reflexx.com)