


**INFORMATIVE NOTE**
**Name: REFLEXX N85**
**Date: 12/04/2024 Rev: 06**
**MANUFACTURER**

 REFLEXX SpA  
 Via Passeri, 2 46019 VIADANA (MN)  
 Tel.: + 39 0375 758891 Fax + 39 0375  
 464504 [info@reflexx.com](mailto:info@reflexx.com)  
[www.reflexx.com](http://www.reflexx.com)

**Pack**

Cod./Size	Box	Carton
N85/M	8032891633793	8032891638798
N85/L	8032891633830	8032891638835
N85/XL	8032891633816	8032891638811
N85/XXL	8032891633854	8032891638859
N85B/M	8032891634059	8032891639054
N85B/L	8032891634066	8032891639061
N85B/XL	8032891634073	8032891639078
N85B/XXL	8032891634080	8032891639085

**DESCRIPTION** ■ Disposable glove powder-free, not sterile ■ **Material:** nitrile ■ **Shape:** ambidextrous; five fingers, beaded cuff ■ **Colour:** orange or black ■ **Feature:** externally “goose bump” textured ■ **Packing:** 50 pieces/box; 10 boxes/carton ■ ISO 2859: AQL 1.5 (level G1)

**CONFORMITY AND PERFORMANCES**
**Personal Protective Equipment cat III**

0465 ANCI\_CIMAC, Via Aguzzafame, 60/b 27029 Vigevano (PV)

Reg EU 2016/425 - EN ISO 374-1:2016 TYPE B - EN ISO 374-5:2016 - EN 374-2:2014 - EN 420:2003 + A1 2009



Reg CE 1935/2004 – Reg UE 10/2011 - DM 21/03/1973

Reg CE 1935/2004 – Reg EU 10/2011 - DM 21/03/19. To such use carefully read the declaration of suitability for food contact (DICA).

 \* **declared only for the orange type**
**DIMENSIONS**

Cod./Size	Length mm	Width mm
N85/ M (8) N85B/ M (8)	250	90-99
N85/ L (9) N85B/ L (9)	260	100-109
N85/ XL (10) N85B/ XL (10)	270	110-119
N85/ XXL (11) N85B/ XXL (11)	280	120-129


Thickness mm	Palm	0.15-0.20
--------------	------	-----------

EN 455-2:2009/A1	Tensile strength	≥6 N
------------------	------------------	------

**EN ISO 374-1:2016**

Chemical compound:	COD	EN 16523-1:2015 TYPE B		EN 374-4:2013
		Performance level	Breakthrough time	% degradation
Sodium hydroxide 40%	K	2	>30'	11.9%
Hydrogen peroxide 30%	P	3	>60'	13.7%
Formaldehyde 37%	T	3	>60'	2.9%
n-Heptane	J	2	>30'	22.2%

EN ISO 374-1:2016/Type B



EN ISO 374-5: 2016



RISKS FROM WHICH THE DPI IS DESIGNED TO PROTECT:

COD		Protection
K Sodium hydroxide 40%	CAS:1310-73-2	Protection from use of inorganic bases
P Hydrogen peroxide 30%	CAS: 7722-84-1	Protection from use of peroxides
T Formaldehyde 37%	CAS: 50-00-0	Protection from the use of aldehydes
J n-Heptane	CAS: 142-82-5	Protection from the use of saturated hydrocarbons

**RECOMENDATIONS OF USE** Multi-purpose glove to protect the hand and/or the handled product, can also be used for general uses ■ The glove is suitable for contact with pure chemical products for a period of less than the permeation time indicated in the chemical resistance table ■ Results obtained under laboratory conditions may not reflect actual working time

**WARNINGS** Before using the product, try it in the real condition of use ■ The glove is a device subject to wear, so frequent replacement is recommended and in all cases where imperfections are observed ■ Wear gloves with dry, clean hands ■ In rare cases where transient hypersensitive skin reactions may occur, discontinue use and consult a doctor

**STORAGE AND DISPOSAL** Store the gloves in their original packaging in a dry place. ■ Do not expose to direct sunlight and heat sources. ■ The expiry date is 3 years from the date of manufacture indicated on the sales package, and is valid for the product properly stored ■ Gloves must be disposed of according to local regulations ■ Declaration of Conformity on [www.reflexx.com](http://www.reflexx.com)